

Product: DR-102 – Optical refractometer

1. Identification of the product

Product: Optical refractometer
 Catalog No.: DR-102
 28,0 62,0 %Brix



The refractive index, of a liquid containing sugar or other solids, and' proportional to its concentration. Refractometers manuals are based on this principle and constructed to be quickly and easily used to measure the proportion of sugar, present in the liquid.

Are therefore various fields of use:

In industries for quality control of the raw materials used and of the final product for fruit, honey, drinks, juices etc.

For wine producers, wineries, analytical laboratories.

The refractometer **DR-102**, used for the solutions that are high in sugar as concentrated fruit juice, jam, ketchup and mix d 'egg.

The refractometer **DR-102** is supplied ready to 'use with calibration kit and padded case.

Technical data:

Parameter	Scale	Resolution	Accuracy
% Sugar	28-62%	0,2%	± 0,2

Note: We reserve the right to change, improve and change the products shown



Keep out of reach children



Dispose of by using the appropriate waste disposal authorities.
 Do not dispose in environment